Reflect your health

**No blood, no pain, no risk, no delay –
haemospect and bilispect**

INNOVATIVE NON-INVASIVE QUANTITATIVE

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MBR Optical Systems is the world's first manufacturer of non-invasive blood meters, which use reflection spectroscopy to quantitatively measure hemoglobin and bilirubin levels. With MBR meters, measurements are fast, reliable and painless.

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A sensor head placed on the skin projects a white light into the underlying tissue via a waveguide. Some of the projected light is absorbed by the various components of tissue, while some of it is reflected. Another waveguide transmits the light reflected as a result of the physical conditions back to the device. A spectrometer breaks the light down into its separate wavelengths and an electronic evaluation unit analyses it. The resulting data is then processed using an algorithm developed by MBR GmbH & Co. KG and visualised on the display of the device. Individual measurements and continuous measurements can also be carried out. The device is powered by batteries. A matching charger is included. With fully charged batteries, the device can measure for at least five hours in continuous operation.

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haemospect® system is a highly sensitive measurement device, designed for non-invasive measurement of blood hemoglobin levels.

The **haemospect®** is highly suitable for use as a fast and cost-effective means of blood donor preselection.

Key benefits:

Blood Donor Centre

- Fast, cost-effective donor selection
- Rapid, non-invasive, painless measurement of blood hemoglobin levels with no risk of infection

Anaesthesia

- Monitoring of Hbg levels during the surgery and on ICU
- Enables better use of blood products

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bilispect® system is a highly sensitive measurement device designed for non-invasive measurement of bilirubin levels.

The fast, safe and careful measurement of the bilirubin value in newborns is essential for an optimal therapy. Unlike other non-invasive bilirubin measurement devices **bilispect®** measures the bilirubin in the skin as well as the bilirubin in the blood.

Key benefits:

- reliable results to inform correct therapy and ongoing treatment decisions
- long-term monitoring during the phototherapy using a button sensor
- appropriate for all newborns
- without needle, painless, no stress for infants
- no disposable materials
- minimization of operating costs

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Company

MBR Optical Systems GmbH & Co. KG was established in 2007 by three shareholders:

Prof. Dr. med. Michael Schietzel is a radiologist/ oncologist and was head physician of the Witten/Herdecke University Community Hospital. After serving as an expert for the German Federal Ministry of Education and Research for many years, he founded the "Herdecke Cancer Research Association". Since the 1980s, he has dedicated himself to the study of measurement by means of reflection spectroscopic processes in the qualitative and quantitative analysis of substances.

Dr. rer. med. Holger Jungmann worked as a mathematician and physicist in basic research at the Max-Planck-Institute. With Prof. Schietzel, he is an architect of the reflection spectroscopic processes applied here and is responsible for technology development for MBR devices.

Dipl. Ing. Heinz Schmersal is managing partner of the Schmersal Group, which specialises in the development and production of highly sensitive safety engineering systems used in industrial, hygienic and medical applications. He is responsible for the production of the MBR-Products.

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United Kingdom (inkl. Nordirland)

L.In.C. Medical Systems Ltd
Loddington House, Main Street
Loddington Leicestershire LE7 9XE
UK
www.linc-medical.co.uk

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Taiwan

Amazing Medical Inc.
3F., No. 13, Lane 13, Lyn-Yi St.
Taipei 100
Taiwan, R.O.C.
www.ammeditech.com

Reflect your health

Hölker Feld 5
D-42279 Wuppertal
Tel. +49 (0)202/64 74-550
Fax +49 (0)202/64 74-560
mbr@mbr-os.com
www.mbr-optical-systems.com

INNOVATIVE NON-INVASIVE QUANTITATIVE

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Stradis-Med GmbH
MicroGene Diagnostic Systems (P) Ltd.