Reflect your health

haemospect[®]

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Friedrich Nietzsche (1844-1900), dt. Philosoph

World premiere

“Today, this morning’s utopia is this
afternoon’s reality.”

Friedrich Nietzsche (1844-1900), German philosopher

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MBR Optical Systems GmbH & Co. KG has successfully developed a highly innovative instrument for the non-invasive, quantitative measurement of haemoglobin.

The **haemospect®** mobile handheld device is a diagnostic system that tremendously improves the quality of the patient's life – painlessly and without the need to draw blood. This means the right therapy decisions can be taken far faster than with previous systems.

These represent outstanding benefits, particularly for use with neonates and premature neonates, in emergency medicine and with burn victims.


Further innovative applications of the **haemospect®** include the measurement of vasomotion, tissue water content and capillary diameter.

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
Haemoglobin O2-Haemoglobin Vasomotion

Blood is principally a transport medium. Haemoglobin, the red colorant in red blood cells, transports vitally important oxygen from the lungs into all body tissues. Even a short interruption of the oxygen supply can lead to life-threatening damage of organs, or even death. This is why haemoglobin plays such a crucial role in clinical diagnostics. But the adequate loading of haemoglobin with oxygen is not the only mechanism that renders blood necessary for life: the haemoglobin in the red blood cells then also has to make its way into the tissues. This requires a vascular system, with its own characteristic movement (vasomotion), to transport the blood to the vicinity of the organs that can absorb the oxygen. The **haemospect®** measures both O2 haemoglobin and vasomotion, rapidly and painlessly.

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Principle of operation

A sensor head placed on the skin projects a white light into the underlying tissue via a waveguide.

Some of the projected light is absorbed by the various components of tissue, while some of it is reflected.

Another waveguide transmits the light reflected as a result of the physical conditions back to the device.

A spectrometer breaks the light down into its separate wavelengths and an electronic evaluation unit connected to the system analyses it.

The resulting data is then processed using an algorithm developed by MBR and visualised on the display of the device in the form of quantitative values measured for the parameters described above.

Individual measurements and continuous measurements can also be carried out online. The device is powered by batteries.

A matching intelligent charger is included, so the batteries will not be overloaded.

With fully charged batteries, the device can measure for at least five hours in continuous operation.

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Benefits at a glance

- Painless measurement of haemoglobin values without the need to draw blood
- Measurement results in seconds
- Documentation of the measured values in the device
- Hygienic, easy, reliable operation
- Rapid application of results for therapy
- Mobility and versatility through battery-powered operation
- No consumables, no laboratory
- Extremely attractive economic benefits, rapid amortisation of investment cost

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Company

MBR Optical Systems was established in 2007 by three shareholders:

Prof. Dr. med. Michael Schietzel is a radiologist/ oncologist and was head physician of the Witten/Herdecke University Community Hospital. After serving as an expert for the German Federal Ministry of Education and Research for many years, he founded the “Herdecke Cancer Research Association“. Since the 1980s, he has dedicated himself to the study of measurement by means of reflection spectroscopic processes in the qualitative and quantitative analysis of substances.

Dr. rer. med. Holger Jungmann worked for many years as a mathematician and physicist in basic research at the Max-Planck-Institute. Together with Prof. Schietzel, he is an architect of the reflection spectroscopic processes applied here and is responsible for technology development for MBR devices.

Dipl. Ing. Heinz Schmorsal is managing partner of the Schmorsal Group, which employs 1,200 people worldwide and specialises in the development and production of highly sensitive safety engineering systems used in industrial, hygienic and medical applications. MBR benefits from his knowledge in the design of manufacturing, production and quality systems.

Prof. Dr. med.
Michael Schietzel

Dr. rer. med.
Holger Jungmann

Dipl. Ing.
Heinz Schmorsal

Reflect your health

Hölker Feld 5
D-42279 Wuppertal
Tel. +49 (0)202/64 74-550
Fax +49 (0)202/64 74-560
mbr@mbr-os.com
www.mbr-optical-systems.com

INNOVATIVE NON-INVASIVE QUANTITATIVE